

cardiac myocytes with AFM. The contractility of cardiac myocytes was measured by using a contact-mode operation, in which the deflection of cantilever with a colloidal probe that contacted the cell surface was kept a constant value with an electronic feedback loop. A force mapping technique was employed to measure the spatial-dependent of the contractility of spontaneously beating cardiac myocytes. We succeeded to measure the amplitude and the frequency of spontaneously beating cardiac myocytes with the AFM technique, and it was observed that the dynamics of the cells was kept in a steady state under appropriate conditions. Interestingly, the cells distended in a region around the center of cells while they exhibited a contraction in a peripheral region of cells. Such a spatial-dependent contractility was observed as the beating was externally perturbed with chemicals. This work is partially supported by the GCOE Program from MEXT of Japan.

**2047-Pos Board B17  
Mechanical Activity At Focal Adhesion Sites**

**Gerald A. Meininger**, Zhe Sun, Shaoxing Huang, Michael A. Meininger.  
University of Missouri, Columbia, MO, USA.

We investigated whether mechanical force applied to extracellular matrix proteins (ECM)-integrin focal adhesion sites would induce mechanical activity characteristic of specific ECM type. We used atomic force microscopy (AFM) to apply forces to ECM adhesion sites on vascular smooth muscle cells (VSMC) isolated from resistance arterioles. The tip of the AFM probes were fused with a borosilicate bead (5  $\mu\text{m}$  diameter) coated with fibronectin (FN), collagen type-I (CNI), collagen type-IV (CNIV), laminin (LN) or vitronectin (VN). ECM-coated beads induced clustering of  $\alpha 5$  and  $\beta 3$  integrins and actin filaments at sites of bead-cell contact indicative of focal adhesion formation. Step increases of an upward (z-axis) pulling force (800~1600 pN) applied to the bead-cell contact site for FN specific focal adhesions induced a force-generating response from the cells resulting in a downward pull by the cell. Depolymerization of the actin cytoskeleton with cytochalasin D blocked whilst stabilization of the actin cytoskeleton with jasplakinolide enhanced this micromechanical event. Myosin light chain kinase inhibition (ML7) and an inhibitor of cSrc tyrosine kinase (PP2) also blocked the response. Furthermore, inhibitory antibodies to  $\alpha 5$  and  $\beta 3$  integrins blocked the micromechanical cell event in a concentration-dependent manner. Similar experiments with CNI, CNIV, VN, or LN failed to induce micromechanical events. Our results demonstrate that mechanical force applied through FN at single focal adhesion sites induces a micromechanical event that is actin, myosin light chain kinase and tyrosine kinase dependent. Importantly, the data illustrate that there are different mechanical characteristics for focal adhesions formed by different ECM proteins. FN appears of particular relevance in its ability to induce a force-generating reaction from sites of focal adhesion in VSMC in response to applied forces.

**2048-Pos Board B18  
Cells and Gels: A Comparison of Indentation Behavior**

David C. Lin<sup>1</sup>, David I. Shreiber<sup>2</sup>, **Emilios K. Dimitriadis**<sup>1</sup>, Ferenc Horkay<sup>1</sup>.  
<sup>1</sup>National Institutes of Health, Bethesda, MD, USA, <sup>2</sup>Rutgers, The State University of New Jersey, Piscataway, NJ, USA.

Understanding and modeling the mechanical behavior of biological systems requires high quality experimental data to be acquired and analyzed. Conventional indentation tests assume homogeneous linear elastic material properties, and therefore are of limited use in characterizing inhomogeneous biological tissues that generally exhibit complex nonlinear elastic response. Despite their shortcomings, the Hertzian models commonly applied in the analysis of the indentation of soft materials are accurate for determining the quasi-static elastic properties of many synthetic gels when applied within their restrictions, including small strain deformation. The lack of practicable, nonlinear elastic contact mechanics models has compelled the application of the Hertz theory to the indentation of biological soft matter, including cells, often with erratic results. We evaluated the accuracy and limitations of the Hertz equation of spherical contact as it pertains to the AFM indentation of chemically crosslinked poly(vinyl alcohol) gels and murine chondrocytes. We then derived and validated via numerical and experimental methods, nonlinear elastic contact equations based on different hyperelastic strain energy functions. The data were reanalyzed using the new models. The results revealed that the linear elastic limit of the cells is dramatically smaller than that of the synthetic gels. All hyperelastic models considered capably described the mechanical response of the gels while only the Fung and Ogden phenomenological models proved suitable for the chondrocytes, which exhibited pronounced strain stiffening even at small deformations. We propose the use of these mathematically simple alternatives to the Hertz theory for modeling the indentation of intrinsically nonlinear soft materials, including live cells and soft tissues, at strains exceeding the Hertzian regime.

**2049-Pos Board B19**

**Physical Properties of Native Biofilm Cells Explored by Atomic Force Microscopy**

Catherine B. Volle<sup>1</sup>, Megan Ferguson<sup>2</sup>, Kathy Aidala<sup>1</sup>, Eileen M. Spain<sup>3</sup>, **Megan E. Nunez**<sup>1</sup>.

<sup>1</sup>Mount Holyoke College, South Hadley, MA, USA, <sup>2</sup>State University of New York, New Paltz, NY, USA, <sup>3</sup>Occidental College, Los Angeles, CA, USA. Biofilms are complex microbial communities that grow at interfaces. Bacteria in biofilms are phenotypically different than their planktonic (free swimming) relatives; they adapt to the communal, sessile lifestyle by optimizing their motility, adhesion, and metabolism. We used Atomic Force Microscopy (AFM) to directly probe the physical properties of native bacterial cells in simple biofilm communities and demonstrated that widely dissimilar biofilm-colonizing cells all have a high cellular spring constant, indicating that they are quite stiff. However, the lab strain *E. coli* ML35 that does not form robust biofilms is much less stiff, hinting that stiff bacteria may preferentially colonize surfaces in the early stages of biofilm formation. Adhesive forces between the retracting AFM tip and bacterial cells vary between cell types in terms of the force components, the distance components, and the number of adhesion events, reflecting differences in associated extracellular polymeric substances (EPS), pili, and flagella. Because biofilms are dynamic, robust, and challenging to control or destroy, potential removal agents are of great interest in industrial, medical, and agricultural settings. We also explored the changes that occur in *E. coli* biofilm cells as they are devoured by *Bdellovibrio bacteriovorus*, a bacterial predator of other bacteria. Invaded prey cells, called bdelloplasts, undergo substantial chemical and physical changes that we probed directly with the AFM tip. Bdelloplasts are significantly shorter than uninvaded *E. coli* biofilm cells, and prey cells clearly lose elasticity after invasion by *B. bacteriovorus* predators. On average, the spring constant of uninvaded *E. coli* cells was three times stiffer than that of bdelloplasts. The retraction portions of the force curves indicate that bdelloplasts adhere to the AFM tip with larger pull-off forces as compared to uninvaded *E. coli*. Thus, AFM provides provocative information about native biofilms.

**2050-Pos Board B20**

**Viscoelastic Indentation of Extremely Soft Biological Samples**

**William J. Adams**<sup>1</sup>, Megan L. O'Grady<sup>1</sup>, Kaustabh Ghosh<sup>2</sup>, Ashley D. Gibbs<sup>2</sup>, Nicholas A. Geisse<sup>1</sup>, Donald E. Ingber<sup>2,1</sup>, Kevin Kit Parker<sup>1</sup>.  
<sup>1</sup>Harvard School of Engineering and Applied Sciences, Cambridge, MA, USA, <sup>2</sup>Children's Hospital and Harvard Medical School, Boston, MA, USA.

Though tissue mechanics serve an important biological function, as in embryonic organ development, a paucity of experimental methods exists to measure these interactions. We used an atomic force microscope (AFM) to measure local mechanical properties of unfixed cryosections of mouse embryonic tooth, which served as a model for organ morphogenesis. AFM is commonly used for tissue indentation; however, most of these studies analyze adult tissues that are mechanically more robust than embryonic tissues, and thus existing experimental protocols do not address the technological limitations of AFM indentation into extremely soft materials (<10 Pa). Importantly, in the range of small applied forces (100 pN), environmental noise, reaction time and drift of the AFM cantilever can apply small, but significant, forces to the sample that affect analyses. Therefore, these artifacts should be incorporated into the indentation load history for viscoelastic analysis. To measure the viscoelastic creep function of embryonic tooth tissue, we performed load-prescribed indentations with an AFM that produced time-dependent load-indentation curves. Our analyses revealed dramatic deviations (>50%) between the software-prescribed and actual applied loads that arise from tissue deformation and the artifacts just mentioned. Since there is no closed-form solution to the integral of the noisy load history, analysis of the indentation data is more complicated. To address this, we took a numerical approach, which produced an average correction of 11% in the instantaneous shear modulus, 10% in the infinite-time shear modulus and 19% in the creep function time constant of relaxation in data from 20 indentations. In addition, this approach was able to successfully resolve subtle differences (<10 Pa) in local tissue stiffness. Thus, our method produces a more sensitive and accurate viscoelastic measurement of extremely soft embryonic tissues, which may greatly aid in uncovering the micromechanical determinants of embryonic tissue development.

**2051-Pos Board B21**

**Statistics of Cell Rheology Measured by AFM**

Shinichiro Hiratsuka, Yusuke Mizutani, Koichi Kawahara, Hiroshi Tokumoto, **Takaharu Okajima**.  
Hokkaido Univ, Sapporo, Japan.

The atomic force microscopy (AFM) proved to be a useful method for measuring the rheological behaviors of living cells. The force-modulation mode of AFM allowed us to measure the complex shear modulus,  $G^*(\omega)$ , quantitatively in a wide frequency range. One of the advantages of the force modulation is that